[ENGLISH](#)

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Sydney, NSW 2000

Australia

EUROPEAN UNION SPECIFIC INFORMATION

Tacx B.V.

De Boeg 2

2343 HK Oegstgeest The Netherlands

Contact: medicalcompliance@garmin.com

SWITZERLAND SPECIFIC INFORMATION

MedEnvoy Switzerland

Gotthardstrasse 28

6302 Zug

Switzerland

Users must read instructions before use.

1. INDICATIONS FOR USE (U.S., U.S. Territories, Australia, Canada, and New Zealand)

The Garmin ECG app is a software-only, mobile medical application intended for use with compatible Garmin smartwatches to create, record, store, transfer, and display a single-channel electrocardiograph similar to a Lead I ECG. The ECG app determines the presence of atrial fibrillation (AFib) or sinus rhythm (SR) on a classifiable waveform. The ECG app is not recommended for users with other known arrhythmias.

The ECG app is intended for over-the-counter (OTC) use. The ECG data displayed by the ECG app is intended for informational use only. The user is not intended to interpret or take clinical action based on the device output without consultation of a qualified healthcare professional. The ECG waveform is meant to supplement rhythm classification for the purposes of discriminating AFib from normal SR, and it is not intended to replace traditional methods of diagnosis or treatment. The ECG app is not intended for use by people under 22 years old.

2. INTENDED PURPOSE (European Union, Switzerland, and United Kingdom)

The Garmin ECG app is a software-only, mobile medical application intended for use with compatible Garmin smartwatches to create, record, store, transfer, and display a single-channel electrocardiograph similar to a Lead I ECG. The ECG app determines the presence of atrial fibrillation (AFib) or sinus rhythm (SR) on a classifiable waveform. The ECG app is not recommended for users with other known arrhythmias. The ECG app is not intended for use by people under 22 years old.

The ECG app is intended for over-the-counter (OTC) use. The user is not intended to interpret or take clinical action based on the device output without consultation of a qualified healthcare professional. The ECG waveform is meant to supplement rhythm classification for the purposes of discriminating AFib from normal SR, and it is not intended to replace traditional methods of diagnosis or treatment.

Target Population and Intended Users

The ECG app is intended for home use by people aged 22 years and over. There is no specific clinical condition required to use the ECG app.

3. CONTRAINDICATIONS, GENERAL WARNINGS, AND PRECAUTIONS

Please observe the following CONTRAINDICATIONS.

Use of the ECG app in these situations could cause harm.

DO NOT use if you are younger than 22 years old.

DO NOT use if you have a known arrhythmia other than AFib.

DO NOT use if you have a cardiac pacemaker, implanted cardiac defibrillator, or other implanted electronic devices.

DO NOT use for clinical diagnosis of heart-related conditions. ECG app results are NOT a diagnosis.

DO NOT use for continuous, real-time monitoring of heart rhythm.

DO NOT use if you are pregnant. The ECG app has not been tested on individuals who are pregnant.

Please observe the following WARNINGS.

Failure to heed these warnings could cause harm.

The ECG app CANNOT check for signs of a heart attack. If you believe you are having a medical emergency, call emergency services.

The ECG app CANNOT detect blood clots, stroke, or other heart conditions, including high blood pressure, congestive heart failure, high cholesterol, and other arrhythmias.

DO NOT change your medication or dosage without first talking to your doctor. ECG app results may not be accurate for people taking certain medications or substances that affect heart rate or blood flow.

DO NOT use in place of any treatment prescribed by your doctor.

ECG app results are only potential findings. DO NOT interpret or take clinical action based on ECG app results without first consulting a qualified healthcare professional.

DO NOT use with a cardiac pacemaker, implanted cardiac defibrillator, or other implanted electronic devices.

A sinus rhythm result only applies to that ECG recording and DOES NOT mean your heart always beats normally. A sinus rhythm result DOES NOT guarantee you are not experiencing an arrhythmia or other health conditions. You should notify your doctor if you detect possible changes in your health.

Even if you have been diagnosed with AFib, your heart may not always beat irregularly. Detecting AFib can be difficult, even in a doctor's office.

Please follow the PRECAUTIONS listed below.

The ECG app may not be able to provide accurate results if you do not follow these precautions.

DO NOT use the ECG app during physical activity.

DO NOT use the ECG app when you are in close vicinity to strong electromagnetic fields (examples: electromagnetic anti-theft systems, metal detectors).

DO NOT use the ECG app during medical procedures, including magnetic resonance imaging (MRI).

DO NOT use the ECG app when your Garmin watch is outside of the temperature range of 50 °F to 113 °F (10 °C to 45 °C) or humidity range of 10% to 95%.

DO NOT use the ECG app while you or the Garmin watch is submerged in water.

Wear your Garmin watch as instructed. The ECG app can only be used on compatible Garmin smartwatches. For information on which Garmin products are compatible with the ECG app, go to garmin.com/ECG.

Make sure your thumb and index finger used to record an ECG are uncovered and free of any bruises, scars, or cuts.

Certain physiological conditions can prevent some people from having a strong enough signal for the ECG app to detect and analyze.

For security reasons, Garmin recommends you maintain control of your Garmin watch and smartphone at all times to prevent unauthorized access to your health data. Garmin also recommends you enable the device lock feature of your Garmin watch. You should enable PIN, password, or fingerprint functionality on your smartphone following the manufacturer's instructions to further prevent unauthorized access. Do not use your smartphone or Garmin watch on unsecured Wi-Fi® networks. You should also secure your Garmin Connect™ account using a strong password and complete setup of multi-factor authentication.

4. ECG APP SETUP

The ECG app can only be used on a compatible Garmin watch that has been set up with your Garmin Connect account. The ECG app requires use of the Garmin Connect app, which is available to download from the app store on your smartphone. For information about required Garmin watch and Garmin Connect app software versions, see the below section **Required Software Versions**.

The ECG app is only available in select regions. For ECG app regional availability and general support, go to garmin.com/ECG.

To use the ECG app on your Garmin watch, you must first complete the setup in the Garmin Connect app.

1. Open the Garmin Connect app on your smartphone.
2. Select the device image at the top of the screen.
Note: If you have more than one device paired with the Garmin Connect app, select the **...** menu, select **Garmin Devices**, and then select the device.
3. Select **Finish Setup > Set Up ECG App**.
4. Follow the on-screen instructions.
5. After the setup is complete, you can record an ECG using the ECG app on your Garmin watch.

5. RECORDING AN ECG ON YOUR GARMIN WATCH

Follow these instructions to record an ECG on your Garmin watch.

1. Wear the watch just above your wrist bones toward your elbow on the wrist you selected during setup.
Note: The watch should be snug but comfortable.
2. Open the ECG app on your watch and follow the on-screen instructions.
3. Sit comfortably, rest your arm and wrist on a table, and keep still.
4. Place your thumb and index finger on the metal ring around the watch face to start a recording.
5. Keep your fingers in place and remain still for the entire 30-second recording.
6. After the recording is complete, you can view the ECG app result on your watch. Please see below for more information about results you might see.

6. VIEWING YOUR ECG RECORDING IN THE GARMIN CONNECT APP (OPTIONAL)

If you chose to upload your health and fitness data during setup, then you can view your ECG recordings in the Garmin Connect app on your smartphone.

1. From the Garmin Connect app, select the **...** menu.
2. Select **Health & Performance (or Health Stats) > ECG**.
3. Choose the ECG recording to view more information.

7. HOW IT WORKS

The ECG app uses sensors on your Garmin watch to record the electrical signals that control how your heart beats. This recording is known as an electrocardiogram, or ECG. The ECG app analyzes the recording to get your heart rate and detect signs of an irregular heart rhythm called AFib. AFib rhythms occur when the upper and lower chambers of your heart are not beating in sync.

8. RESULTS YOU MIGHT SEE

After you complete an ECG recording, you will see one of the following results in the ECG app.

Sinus Rhythm

A sinus rhythm result means your heartbeat pattern appears normal. It indicates that the upper and lower chambers of your heart are beating in sync.

Note: A sinus rhythm result only applies to that specific recording and DOES NOT mean your heart always beats normally. Notify your doctor if you detect any possible changes in your health.

Atrial Fibrillation (AFib)

An AFib rhythm result means your heartbeat pattern appears irregular. AFib is an irregular heart rhythm caused by your heart's upper and lower chambers not beating in sync. You should speak with your doctor if you receive an AFib rhythm result. Left untreated, AFib can lead to serious health consequences. The irregular heart rhythm causes improper blood flow in the heart, which can lead to stroke, heart failure, and/or other medical problems.

AFib can be temporary or persistent, but it is a manageable condition when treated by a doctor and with medication. Many people with AFib live healthy lives. Sometimes people with AFib do not have symptoms. However, others may experience one or more of these common symptoms:

- Rapid heartbeat (pounding, fluttering)
- Skipped heartbeat
- Fatigue
- Shortness of breath
- Chest tightness or pain
- Fainting
- Dizziness

Heart Rate Too High or Too Low

If your heart rate is over 120 beats per minute (bpm) or under 50 bpm, the ECG app cannot check for AFib heart rhythms. In these cases, the result will be inconclusive, and your heart rhythm will not be classified. If you get this result repeatedly or you do not feel well, then you should speak with your doctor.

A heart rate below 50 bpm can be caused by taking some medications or by certain heart conditions. Some people with high fitness levels may also have a heart rate below 50 bpm.

A heart rate above 120 bpm is common while exercising and for a short time afterward. It can also be caused by high stress levels, an infection, dehydration, alcohol use, and certain heart conditions, including AFib.

Inconclusive

An inconclusive result means your heart rhythm could not be classified. This could happen for several reasons, such as moving too much or poor skin contact during the ECG recording. Try resting your arm and wrist on a table and make sure the watch is snug on your wrist. If you get this result repeatedly or if you do not feel well, then you should speak with your doctor.

Here are some steps to help improve your recording:

- Rest your arm and wrist on a table and keep still during the recording.
- Make sure the Garmin watch is snug on your wrist. The back of the watch must be touching your arm.
- Only place your thumb and index finger on the metal ring around the watch face and keep them in place for the full recording.
- Make sure you are wearing the watch on the wrist you chose during setup.
- Make sure both your wrist and the watch are clean and dry.
- Stay clear of other electronics to avoid interference.

Other factors that can cause an inconclusive result:

- A heart rate between 100 and 120 bpm, and the rhythm was not AFib.
- A pacemaker or other implanted cardioverter device (ICD).
- Other heart conditions that the ECG app cannot recognize.
- Not enough electrical signal in the recording.

9. SAFETY AND PERFORMANCE

The ECG app's ability to accurately detect AFib and sinus rhythms in an ECG recording was validated in a clinical study involving approximately 590 subjects. ECG app rhythm classifications were compared to 12-Lead ECG rhythm classifications performed by board-certified cardiologists. The ECG app was able to correctly identify AFib 99.5% of the time and correctly identify sinus rhythm (SR) 100% of the time in the recordings that could be classified.

During this study, the ECG app determined 11.5% of recordings were inconclusive, including those which were deemed of too poor quality for analysis. When including these inconclusive recordings, the probability that the ECG app would return an AFib result for a subject in AFib was 86.5%, and 91.1% for a SR result for subjects in sinus rhythm. Real-world performance may have a higher rate of inconclusive and poor recording results.

Accuracy of the ECG app PDF reports was assessed by comparing to a simultaneously recorded standard Lead I ECG. Key characteristics of the ECG waveforms such as PR and RR intervals, QRS duration, location, and amplitude, and P wave location and amplitude were compared and found to be statistically equivalent within a 1 millimeter non-inferiority margin. Comparison of the board-certified cardiologists' rhythm classifications of simultaneously recorded 12-Lead ECG and ECG app recordings had 96% agreement. No adverse events were observed during the clinical study.

10. CLINICAL BENEFITS

The ECG app provides a convenient means for at-home ECG recording during times of perceived symptoms. This opportunity increases the likelihood of early AFib detection, which allows you to seek medical attention earlier. Early clinical intervention could potentially delay, reduce, or control long-term complications of AFib. Sharing your ECG app results with a healthcare professional may improve the efficiency of your medical evaluation.

11. TROUBLESHOOTING

If you experience problems using the ECG app, refer to this troubleshooting guide.

Note: Any serious incidents which occur in relation to the ECG app should be reported to the manufacturer. For the European Union, serious incidents should also be reported to the European Authorized Representative and the competent authority of the member state in which you reside.

Problem: I can't activate the ECG app on my compatible Garmin watch.

Solution:

- Make sure the ECG app is available for your Garmin watch. Go to garmin.com/ECG.
- Make sure the ECG app is approved for use in your country. Go to garmin.com/ECG.
- Make sure you have location permissions enabled for the Garmin Connect app on your smartphone. Refer to the manufacturer's instructions to enable location permissions.
- Review and complete the setup in the Garmin Connect app on your smartphone (see the above section **ECG App Setup**).

Problem: I can't find the ECG app on my Garmin watch.

Solution:

- Make sure the ECG app is available for your Garmin watch. Go to garmin.com/ECG.
- Make sure the ECG app is approved for use in your country. Go to garmin.com/ECG.
- Complete the setup in the Garmin Connect app on your smartphone (see the above section **ECG App Setup**).
- The ECG app can be found in the Apps or Activities menu on your Garmin watch. For help finding these menus, please refer to your Garmin watch user manual.

Problem: I can't get the ECG app to create a recording.

Solution:

- Complete the setup in the Garmin Connect app on your smartphone (see the above section **ECG App Setup**).
- Water and sweat can cause a poor recording. Clean and dry your wrist and your Garmin watch.
- Dry skin, cold skin, a hairy wrist, wearing the watch on your wrist bones, or too much movement can prevent the ECG app from recording your ECG.
 - Use moisturizing lotion on your hands and wrist.
 - Remove the watch and rub your wrist gently to warm the skin. Put the watch back on and try again.
 - Wear the watch just above your wrist bones toward your elbow.
 - Remain still during the recording.

Problem: The ECG app says my recording was inconclusive due to a poor quality signal, or I see a lot of noise or interference in my recording.

Solution:

- Sit comfortably, relax, and try to remain still.
- Rest your arm and wrist on a table while you take a recording.
- Check the fit of the Garmin watch on your wrist. The band should be snug. The back of your Garmin watch must be touching your wrist. Wear the watch just above your wrist bones toward your elbow.
- Use your thumb and index finger on the metal ring around the watch face.

- Move away from any electronics that are plugged into a power source to avoid electrical interference.
- Water and sweat can cause a poor recording. Clean and dry your wrist and your Garmin watch.

Problem: The ECG app says my recording was inconclusive, and I do not understand why.

Solution:

Several factors can cause this inconclusive result:

- A heart rate between 100 and 120 bpm, and the rhythm was not AFib.
- Other heart conditions that the ECG app cannot recognize.
- A pacemaker or other implanted cardioverter device (ICD).
- Certain physiological conditions that prevent some people from having enough electrical signal for a recording.

Problem: The ECG waveforms appear upside down.

Solution:

- You may not be wearing the watch on the wrist you selected while setting up your Garmin watch. You can either wear the watch on your other arm or change your user profile wrist selection. To change your wrist selection, please refer to **Setting Up Your User Profile** in your Garmin watch user manual and select **Wrist**.

Problem: I can't find my ECG recordings in the Garmin Connect app.

Solution:

- If you chose to upload your health and fitness data during setup, you can view your ECG recordings in the Garmin Connect app on your smartphone and manually share them by creating a PDF.
- From the Garmin Connect app, select the ⋮ menu.
- Select **Health & Performance** (or **Health Stats**) > **ECG**.

Problem: My Garmin watch storage is full.

Solution:

- If your Garmin watch storage is full, you will not be able to take new ECG recordings. To make room for new ECG recordings, free up space on your Garmin watch by deleting unwanted apps, music, or podcasts.

Problem: I want to stop using the ECG app, or I want to delete all my ECG app results from Garmin Connect.

Solution:

- To remove the ECG app from your Garmin watch and delete any unsynced results, please follow the directions in your Garmin watch user manual to **Delete Data and Reset Settings**.
- To remove your ECG app results from Garmin Connect, go to garmin.com/account and select **Manage Your Data**.

12. RESIDUAL RISKS

The ECG app has many measures in place to help ensure safe use and performance, however you should be aware that some level of residual risk associated with use of the ECG app remains.

Residual Risk	Guidance
A poor-quality ECG signal may result in a failure to detect atrial fibrillation, which could lead to a failure to seek treatment.	Follow the instructions above in Recording an ECG on Your Garmin Watch and Troubleshooting .
An incorrect rhythm classification could lead to a failure to seek treatment or unnecessarily seeking treatment.	If you are not feeling well, or if you doubt the ECG app result, speak to your doctor. If you believe you are having a medical emergency, contact emergency services.
A misinterpretation of ECG app results could lead to inappropriate modification of established treatment plans, failure to seek treatment, and/or unnecessarily seeking treatment.	Consult the above section Results You Might See for a detailed explanation of all ECG app results.
Use of the ECG app to detect conditions other than atrial fibrillation and/or when you have other known arrhythmias could lead to a failure to seek treatment.	The ECG app is not recommended for those with other known arrhythmias. Please observe all information in the above section Contraindications, General Warnings, and Precautions .
A cybersecurity exploit could result in: • Unintentional disclosure of ECG app results; • Erroneous or malicious alteration of ECG app results and/or an inability to use the ECG app that leads to either a failure to seek treatment or unnecessarily seeking medical procedures.	Follow all security recommendations listed in the above section Contraindications, General Warnings, and Precautions .

13. EQUIPMENT SYMBOLS

Manufacturer

European Authorized Representative

Authorized Representative for Switzerland

CE Mark

Medical Device

Users must read instructions before use.

Contraindications. Use of the ECG app in these situations could cause harm.

General warnings. Failure to heed these warnings could cause harm.

General precautions. The ECG app may not be able to provide accurate results if precautions are not followed.

14. REQUIRED SOFTWARE VERSIONS

The following table lists the minimum Garmin watch firmware and Garmin Connect app version required to use the most recent ECG app which has been made available. The ECG app version number may be found on the About page of your Garmin watch and the Settings page in Garmin Connect.

Note: The ECG app may not be available for all watches in all regions. For ECG app regional availability and general support, go to garmin.com/ECG.

Garmin Watch		ECG App	Firmware	Garmin Connect	Android™	iOS®
D2™	D2 Mach 1 Pro	1.1.3	14.68	5.1	9	16
Enduro™	Enduro 3	1.1.3	11.28	5.4	9	16
epix™	epix Pro (Gen 2)	1.1.3	14.68	5.1	9	16
	epix Pro (Gen 2) – Porsche® Edition	1.1.3	15.77	5.1	9	16
fēnix®	fēnix 7 Pro	1.1.3	14.68	5.1	9	16
	fēnix 8	1.1.3	11.28	5.4	9	16
Forerunner®	Forerunner 970	1.1.3	5.02	5.11	9	16
quatix®	quatix 7 Pro	1.1.3	14.68	5.1	9	16
	quatix 8	1.1.3	15.31	5.13	9	16
tactix®	tactix 7 – AMOLED Edition	1.1.3	14.68	5.1	9	16
	tactix 8	1.1.3	12.57	5.7	9	16
Venu®	Venu 2 Plus <i>(United States Only)</i>	1.1.3	11.21	5.1	9	16
	Venu 3	1.1.3	7.07	5.1	9	16

© 2025 Garmin Ltd. or its subsidiaries.

All rights reserved. Under the copyright laws, this manual may not be copied, in whole or in part, without the written consent of Garmin. Garmin®, the Garmin logo, fēnix®, Forerunner®, quatix®, tactix®, and Venu® are trademarks of Garmin Ltd. or its subsidiaries, registered in the USA and other countries. Garmin Connect™, epix™, Enduro™, and D2™ are trademarks of Garmin Ltd. or its subsidiaries. These trademarks may not be used without the express permission of Garmin. Wi-Fi® is a registered trademark of Wi-Fi Alliance Corporation. The Porsche wordmark is a registered trademark of Dr. Ing. h.c. F. Porsche AG. Dr. Ing. h.c. F. Porsche AG is a company registered at Stuttgart, Registriergericht HRB 730623. iOS® is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license by Apple Inc. Android™ is a trademark of Google LLC.

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© 2025 Garmin Ltd. eller dess dotterbolag.